# IndelCheck™ CRISPR/TALEN 插入缺失 检测体系

**Cat. No. ICPE-050**

**Cat. No. ICPE-200**

**Cat. No. TENI-050**

**Cat. No. TENI-200**

**Cat. No. TPCR-050**

**Cat. No. TPCR-200**

**Cat. No. IC007**

**Cat. No. IC008**

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
Did you mean CCR5 as a gene symbol?  
Search Gene for [CCR5](#) as a symbol.

**Results: 1 to 20 of 411**      << First < Prev Page 1 of 21 Next > Last >>

Filters activated: Current only. [Clear all](#) to show 417 items.

Name/Gene ID	Description	Location	Aliases	MIM
<input type="checkbox"/> <a href="#">CCR5</a> ID: 1234	chemokine (C-C motif) receptor 5 (gene/pseudogene) [ <i>Homo sapiens</i> (human)]	Chromosome 3, NC_000003.12 (46370142..46376206)	CC-CKR-5, CCCKR5, CCR-5, CD195, CKR-5, CKR5, CMKBR5, IDDM22	601

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Gene ID: 1234, updated on 17-May-2015

**Summary**

**Official Symbol** CCR5 provided by [HGNC](#)

**Official Full Name** chemokine (C-C motif) receptor 5 (gene/pseudogene) provided by [HGNC](#)

**Primary source** [HGNC:HGNC:1606](#)

**See related** [HPRD:03223](#); [MIM:601373](#)

**Gene type** protein coding

**RefSeq status** REVIEWED

**Organism** [Homo sapiens](#)


**Lineage** Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Primates; Haplorhini; Catarrhini; Hominidae; Homo

**Also known as** CKR5; CCR-5; CD195; CKR-5; CCCR5; CMKBR5; IDDM22; CC-CKR-5

**Summary** This gene encodes a member of the beta chemokine receptor family, which is predicted to be a seven transmembrane protein similar to G protein-coupled receptors. This protein is expressed

**Table of contents**

- Summary
- Genomic context
- Genomic regions, transcripts, and products
- Bibliography
- Phenotypes
- Variation**
- HIV-1 interactions
- Pathways from BioSystems
- Interactions
- General gene information
  - Markers, Clone Names, Homology, Gene Ontology
- General protein information
- NCBI Reference Sequences (RefSeq)

**Variation**

[See variants in ClinVar](#)


[See studies and variants in dbVar](#)

[See Variation Viewer \(GRCh37.p13\)](#)

Genotypes

[See SNP Geneview Report](#)

[See 1000 Genomes Browser \(GRCh37.p13\)](#)

**SNP linked to Gene (geneID:1234) Via Contig Annotation**

The SNP GeneView page only reports human variation on GRCh38. A new [Variation Viewer](#) is available to view the gene CCR5 variations in [GRCh37p13](#) or [GRCh38](#), and will replace SNP GeneView later this year. Please visit the [Help Page](#) or [YouTube](#) for available features and send your comments and suggestions to NCBI [helpdesk](#).

Send  on all gene models to Batch Query  all rs# to file.

**Gene Model (mRNA alignment) information from genome sequence** ↑

Total gene model (contig mRNA transcript): 2

mrna	transcript	protein	mrna orientation	Contig	Contig Label	List SNP
<a href="#">NM_000579.3</a>	plus strand	<a href="#">NP_000570.1</a>	forward	<a href="#">NT_022517.19</a>	GRCh38.p2	<- currently shown
<a href="#">NM_001100168.1</a>	plus strand	<a href="#">NP_001093638.1</a>	forward	<a href="#">NT_022517.19</a>	GRCh38.p2	<a href="#">View snp on GeneModel</a>

Clinical Source  in gene region  cSNP  has frequency  double hit

gene model	Contig Label	Contig	mrna	protein	mrna orientation	transcript	snp count
(contig mRNA transcript):	GRCh38.p2	<a href="#">NT_022517.19</a>	<a href="#">NM_000579.3</a>	<a href="#">NP_000570.1</a>	forward	plus strand	259, coding

Region	Chr. position	mRNA pos	dbSNP rs# cluster id	Heterozygosity	Validation	MAF	Allele origin	3D	Linkout	Function	dbSNP allele	Protein residue	Codon pos	Amino acid pos	PubMed
	<a href="#">46372912</a>	<a href="#">367</a>	<a href="#">rs748244565</a>	N.D.						missense	G	Glu [E]	1	4	
										contig reference	C	Gln [Q]	1	<a href="#">4</a>	
	<a href="#">46372913</a>	<a href="#">368</a>	<a href="#">rs745912425</a>	0.000						frame shift	-	Arg [R]	2	4	
										contig reference	AA	Gln [Q]	2	<a href="#">4</a>	
	<a href="#">46372915</a>	<a href="#">370</a>	<a href="#">rs763192695</a>	0.000						missense	A	Met [M]	1	5	
										contig reference	G	Val [V]	1	<a href="#">5</a>	
	<a href="#">46372916</a>	<a href="#">371</a>	<a href="#">rs766432600</a>	0.000						missense	G	Gly [G]	2	5	
										contig reference	T	Val [V]	2	<a href="#">5</a>	
	<a href="#">46372934</a>	<a href="#">389</a>	<a href="#">rs751603911</a>	0.000						missense	G	Gly [G]	2	11	
										contig reference	A	Asp [D]	2	<a href="#">11</a>	
	<a href="#">46372936</a>	<a href="#">391</a>	<a href="#">rs755442066</a>	0.000						missense	G	Val [V]	1	12	
										contig reference	A	Ile [I]	1	<a href="#">12</a>	
	<a href="#">46372937</a>	<a href="#">392</a>	<a href="#">rs781613165</a>	0.000						missense	C	Thr [T]	2	12	
										contig reference	T	Ile [I]	2	<a href="#">12</a>	
	<a href="#">46372939</a>	<a href="#">394</a>	<a href="#">rs753095965</a>	0.000						missense	G	Asp [D]	1	13	

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Following terms and conditions apply to use of the IndelCheck™ CRISPR/TALEN insertion or deletion detection system. If the terms and conditions are not acceptable, the Product in its entirety must be returned to GeneCopoeia within 5 calendar days. A limited End-User license is granted to the purchaser of the Product. The Product shall be used by the purchaser for internal research purposes only. The Product is expressly not designed, intended, or warranted for use in humans or for therapeutic or diagnostic use. The Product must not be resold, repackaged or modified for resale, or used to manufacture commercial products or deliver information obtained in service without prior written consent from GeneCopoeia. Use of any part of the Product constitutes acceptance of the above terms.

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GeneCopoeia warrants that the Product meets the specifications described in the accompanying Product Datasheet. If it is proven to the satisfaction of GeneCopoeia that the Product fails to meet these specifications, GeneCopoeia will replace the Product. In the event a replacement cannot be provided, GeneCopoeia will provide the purchaser with a refund. This limited warranty shall not extend to anyone other than the original purchaser of the Product. Notice of nonconforming products must be made to GeneCopoeia within 30 days of receipt of the Product. GeneCopoeia's liability is expressly limited to replacement of Product or a refund limited to the actual purchase price. GeneCopoeia's liability does not extend to any damages arising from use or improper use of the Product, or losses associated with the use of additional materials or reagents. This limited warranty is the sole and exclusive warranty. GeneCopoeia does not provide any other warranties of any kind, expressed or implied, including the merchantability or fitness of the Product for a particular purpose.

GeneCopoeia is committed to providing our customers with high-quality products. If you should have any questions or concerns about any GeneCopoeia products, please contact us at 301-762-0888.

© 2015 GeneCopoeia, Inc.