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Reverse Transcription Buffer (10 $\times$ )	2.0 $\mu$ L
25 mM dNTP	1.0 $\mu$ L
RNase Inhibitor	1.0 $\mu$ L
Reverse Transcription Enzyme	1.0 $\mu$ L
Total	20.0 $\mu$ L

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Cycle	Steps	Temperature	Duration	Read
1	Denaturation	95 °C	10 min	off
	Denaturation	95°C	10 sec	off
40	Annealing	60 °C	20 sec	off
	Extension	72 °C	15 sec	on

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Temperature	Interval	Duration	Read
72-95 °C	0.5°C	6 sec/each	on
30 °C		30 sec	off

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## Limited Use License and Warranty

### Limited use license

The following terms and conditions apply to use of the Lenti-Pac™ qRT-PCR Titration Kit (the Product). If the terms and conditions are not acceptable, the Product in its entirety must be returned to GeneCopoeia within 5 calendar days. A limited End-User license is granted to the purchaser of the Product. The Product shall be used by the purchaser for internal research purposes only. The Product is expressly not designed, intended, or warranted for use in humans or for therapeutic or diagnostic use. The Product must not be resold, repackaged or modified for resale, or used to manufacture commercial products or deliver information obtained in service without prior written consent from GeneCopoeia. Use of any part of the Product constitutes acceptance of the above terms.

### Limited warranty

GeneCopoeia warrants that the Product meets the specifications described in the accompanying Product Datasheet. If it is proven to the satisfaction of GeneCopoeia that the Product fails to meet these specifications, GeneCopoeia will replace the Product. In the event a replacement cannot be provided, GeneCopoeia will provide the purchaser with a refund. This limited warranty shall not extend to anyone other than the original purchaser of the Product. Notice of nonconforming products must be made to GeneCopoeia within 30 days of receipt of the Product. GeneCopoeia's liability is expressly limited to replacement of Product or a refund limited to the actual purchase price. GeneCopoeia's liability does not extend to any damages arising from use or improper use of the Product, or losses associated with the use of additional materials or reagents. This limited warranty is the sole and exclusive warranty. GeneCopoeia does not provide any other warranties of any kind, expressed or implied, including the merchantability or fitness of the Product for a particular purpose.

GeneCopoeia is committed to providing our customers with high-quality products. If you should have any questions or concerns about any GeneCopoeia products, please contact us at 020-32068595.