

Transferring an ORF from an EZShuttle[™] Gateway[®] PLUS ORF clone or ORFeome Shuttle Clone to a Destination Vector

Reagents

EZShuttle[™] Gateway PLUS ORF clone or ORFeome Gateway entry clone 5X Reaction Buffer (GeneCopoeia Cat. No. RCBM-1001-020) EZRecombinase[™] LR mix (GeneCopoeia Cat. No. RCBM-1001-020) Linearized destination vector 50 mM TE buffer, pH 8.0 Proteinase K GCI-5α chemically competent *E. coli* cells (GeneCopoeia Cat. No. CC001) S.O.C. Medium LB plates with 100 µg/ml ampicillin

Equipment

Autoclaved assorted pipet tips Pipetman pipets: P2, P20, P200, P1000 Autoclaved 1.5-ml microtubes 17 x 100 mm polypropylene tubes (Falcon 2059) 37°C and 42°C incubators

Procedures

1. LR Reaction

Use the following procedure to perform an LR recombination reaction. For a positive control, use 100 ng (2 µl) of pENTR[™]-eGFP.

1.1 Add the following components to a 1.5-ml microcentrifuge tube at room temperature and mix:

Component	Volume
Entry clone (100–300 ng)	1−10 µl
Destination vector (150 ng/µl)	2 µl
5X Reaction Buffer	4 µl
TE buffer, pH 8.0	to 16 µl

- 1.2 Remove EZRecombinase[™] LR mix from −80°C and thaw on ice for about 2 minutes. Vortex briefly twice (2 seconds each time).
- 1.3 To each sample (Step 1.1 above), add 4 µl of EZRecombinase[™] LR mix to the reaction and mix well by vortexing briefly twice. Microcentrifuge briefly.
- 1.4 Return EZRecombinase[™] LR mix to −80°C storage immediately after use.
- 1.5 Incubate reactions at 25°C for 60 minutes.

www.genecopoeia.com

2. Transformation

2.1 Transform 1 µl of each LR reaction into 50 µl of GCI-5α chemically competent *E. coli* cells. Incubate on ice for 30 minutes. Heat shock cells by incubating at 42°C for 30 seconds. Add 450 µl of S.O.C. Medium and incubate at 37°C for 1 hour with shaking. Plate 20 µl and 100 µl of each transformation onto selective plates.

Note: Any competent cells with a transformation efficiency of >1.0 × 10^8 transformants/µg may be used.

2.2 Transform 1 μl of pUC19 DNA (10 ng/ml) or other suitable plasmid for use as a control into 50 μl of GCI-5α chemically competent *E. coli* cells as described above. Plate 20 μl and 100 μl on LB plates containing 100 μg/ml ampicillin

Expected results

An efficient LR recombination reaction will produce > 5000 colonies if the entire transformation is plated.

LB plates (per liter)

- 10.0 g Bacto-Tryptone
- 5.0 g Bacto-Yeast extract
- 5.0 g NaCl
- 15.0 g Agar

Adjust the pH to 7.0 with NaOH (~200 µl 5M NaOH). Autoclave (keep the top loosened to allow steam to vent) and cool to 50°C before adding antibiotics. Mix well and pour LB media into plates.

The purchase of ORFEXPRESS CLONES IS ACCOMPANIED BY A LIMITED LABEL LICENSE UNDER U.S. PATENTS 5,888,732; 6,143,557; 6,171,861; 6,270,969; 6,277,608 OWNED BY INVIROGEN, CORP. ONLY TO USE THE PRODUCT FOR THE INTERNAL RESEARCH PURPOSES OF THE PURCHASER. FOR INFORMATION ON PURCHASING A LICENSE TO USE THE PURCHASED PRODUCT FOR PURPOSES OTHER THAN THE INTERNAL RESEARCH OF THE PURCHASER, CONTACT THE DIRECT OF LICENSING AT INVIROGEN, CORP., 1600 FARADAY AVENUE, CARLSBAD, CA 92008.

For more information, please contact us:

GeneCopoeia, Inc. 9620 Medical Center Drive, #101 Rockville, MD 20850 Tel: 301-762-0888 Fax: 301-762-3888 support@genecopoeia.com www.genecopoeia.com