Recombinant Human ATG4A Protein Data Sheet

Catalog #	hRP-T3476-EF011
Size	25 μg
Protein Name	Human cysteine protease
Protein Symbol	ATG4A
Original Source	Homo sapiens
Expression System	E.coli
GenBank Accession #	NM 052936.1
Uniprot Accession #	Q8WYN0
Description	Cysteine protease ATG4A belongs to a member of the autophagin protein family. This protein is also designated as a member of the C-54 family of cysteine proteases. Autophagy is the process by which endogenous proteins and damaged organelles are destroyed intracellularly. Autophagy is postulated to be essential for cell homeostasis and cell remodeling during differentiation, metamorphosis, non-apoptotic cell death, and aging. Reduced levels of autophagy have been described in some malignant tumors, and a role for autophagy in controlling the unregulated cell growth linked to cancer has been proposed. Cysteine protease ATG4A cleaves the C-terminal part of either MAP1LC3, GABARAPL2 or GABARAP, allowing the liberation of form I. A subpopulation of form I is subsequently converted to a smaller form (form II). Form II, with a revealed C-terminal glycine, is considered to be the phosphatidylethanolamine (PE)-conjugated form, and has the capacity for the binding to autophagosomes. Preferred substrate is GABARAPL2 followed by MAP1LC3A and GABARAP.
Application	WB, ELISA, IP, antibody production, protein array
Fusion tag	N-His
Peptide Length	413aa(including fusion tag)
Molecular Weight	47.3kDa(including fusion tag)
pI	5.1
Activity	NA
5	Storage buffer: 20mM Tris.Cl, 100mM NaCl, 0.4M Urea, 50% Glycerol, pH8.0. Store at -80°C
Storage	and avoid repeated freeze-thaw cycles.
	67 54 45 35 28
Reference:	 A high-throughput FRET-based assay for determination of Atg4 activity Kinetics comparisons of mammalian Atg4 homologues indicate selective preferences toward diverse Atg8 substrates Reactive oxygen species are essential for autophagy and specifically regulate the activity of Atg4



GeneCopoeia Inc.

9620 Medical Center Drive, Suite 101 Rockville, MD 20850 USA

Phone: 301-762-0888 Toll free: 1-866-360-9531 Fax: 301-762-3888

Web: www.genecopoeia.com

Inquiry: inquiry@genecopoeia.com

 $Technical \ Support: \underline{support@genecopoeia.com}$