


EZShuttle™ Recombination Cloning Expression Vector pEZ-LV201

10258 bp, circular



Vector Features

Name	Start	End
attR1	2	126
ccdB	1224	1529
attR2	1570	1693
SV40 promoter	1701	2068
eGFP	2103	2822
IRES	2836	3412
Puromycin resistance gene	3445	4044
WPE	4058	4649
3'-LTR (delta U3/U5)	4852	5088
pUC Ori (complementary strand)	5547	6217
Ampicillin resistance gene (complementary strand)	6362	7222
5'-LTR (RSV/U5)	7637	8048
RRE	8719	8952
cPPT	9444	9567
CMV Promoter	9664	10258

Sequencing Primers

Primer	Sequence
Forward	5'- gcggtaggcgtgtacggt -3'
Reverse	5'- ctggaatagctcagaggc -3'

GeneCopoeia, Inc.
9620 Medical Center Drive, #101
Rockville, MD 20850
Tel: 301-762-0888 Fax: 301-762-3888
Email: inquiry@genecopoeia.com Web: www.genecopoeia.com

GeneCopoeia products are for research use only

© 2015, GeneCopoeia, Inc