


## EZShuttle™ Recombination Cloning Expression Vector pEZ-P002

11551 bp, circular



### Vector Features

Name	Position
T-DNA right border	54 78
attR1	328 452
lac UV5 promoter	489 519
CmR	573 1232
ccdB	1574 1879
attR2	1920 2044
eGFP	2049 2768
NOS terminator	2777 3186
NOS terminator (complementary strand)	3220 3479
NptII (complementary strand)	3869 4663
NOS promoter (complementary strand)	4676 4982
T-DNA left border	5073 5097

### Sequencing Primers

Primer	Sequence
Forward	5'-CAGGGTTTTCCCAGTCACGAC -3'
Reverse	5'-GAACAGCTCCTCGCCCTTGCTCAC -3'

GeneCopoeia, Inc.  
9620 Medical Center Drive, #101  
Rockville, MD 20850  
Tel: 301-762-0888 Fax: 301-762-3888  
Email: [inquiry@genecopoeia.com](mailto:inquiry@genecopoeia.com) Web: [www.genecopoeia.com](http://www.genecopoeia.com)

GeneCopoeia products are for research use only

© 2015, GeneCopoeia, Inc