


EZShuttle™ Recombination Cloning Expression Vector pEZ-P006

11165bp, circular



Vector Features

Name	Position
T-DNA right border	54 78
attR1	338 462
lac UV5 promoter	499 529
CmR	583 1242
ccdB	1584 1889
attR2	1930 2054
NOS terminator	2070 2479
NOS terminator (complementary strand)	2524 2789
hpt (complementary strand)	2925 3950
NOS promoter (complementary strand)	3978 4264
T-DNA left border	4687 4711

Sequencing Primers

Primer	Sequence
Forward	5'-CAGGGTTTTCCAGTCACGAC -3'
Reverse	5'-GGGAAATTCGAGCTCTAAGC -3'

GeneCopoeia, Inc.
9620 Medical Center Drive, #101
Rockville, MD 20850
Tel: 301-762-0888 Fax: 301-762-3888
Email: inquiry@genecopoeia.com Web: www.genecopoeia.com

GeneCopoeia products are for research use only

© 2015, GeneCopoeia, Inc