


## EZShuttle™ Recombination Cloning Expression Vector pEZ-M09

7087bp, circular



### Vector Features

Name	Position
f1 ori	21 449
SV40 pomoter	476 785
Neomycin resistance gene	859 1653
SV40 poly A	1829 1959
pUC ori(complementary strand)	2401 2989
Ampicillin resistance gene(complementary strand)	3160 4020
CMV promoter	4388 4975
T7 promoter	5019 5037
attR1	5070 5194
lac promoter	5219 5249
CmR	5303 5962
ccdB	6304 6609
attR2	6648 6772
Myc tag	6792 6821
bGH poly A	6835 7062

### Sequencing Primers

Primer	Sequence
Forward	5'-gcggtaggcgtgtacggt-3'
Reverse	5'-gtggcacctccagggtc-3'

GeneCopoeia, Inc.  
9620 Medical Center Drive, #101  
Rockville, MD 20850  
Tel: 301-762-0888 Fax: 301-762-3888  
Email: [inquiry@genecopoeia.com](mailto:inquiry@genecopoeia.com) Web: [www.genecopoeia.com](http://www.genecopoeia.com)

GeneCopoeia products are for research use only

© 2016, GeneCopoeia, Inc