


## pEZX-FR01

6962bp, circular



**MCS1:** SpeI, MluI, SmaI, Sall, PmeI, XbaI, BamHI

### Vector Features

Name	Position
TPS	489 649
MCS1	650 695
hLuc	702 2354
Synthesis polyA	2361 2414
SV40 promoter	2434 2761
SV40 ori	2612 2747
Rluc	2824 3759
IRES	3775 4348
Puromycin resistance gene	4381 4980
SV40 polyA	5071 5292
Kanamycin resistance gene	5413 6222
pUC ori	6286 6959

### Sequencing Primers

Primer	Sequence
Forward	5'-agttacttaagctcgggcc-3'
Reverse	5'-atgtgggcatcgggtgaag-3'

GeneCopoeia, Inc.  
9620 Medical Center Drive, #101  
Rockville, MD 20850  
Tel: 301-762-0888 Fax: 301-762-3888  
Email: [inquiry@genecopoeia.com](mailto:inquiry@genecopoeia.com) Web: [www.genecopoeia.com](http://www.genecopoeia.com)

GeneCopoeia products are for research use only

© 2016, GeneCopoeia, Inc