


pEZX-FR02

6300bp, circular



MCS2: AsiSI, EcoRI, BsiWI, PmeI, XhoI, SpeI

Vector Features

Name	Position
SV40 polyA	6 227
CMV promoter	248 841
Synthetic Renilla Luciferase	867 1802
SV40 polyA	1833 2054
Ampicillin resistance gene	2201 3061
pUC ori	3232 3820
SV40 promoter	3949 4306
chimeric intron	4376 4508
hLuc	4583 6235
MCS2	6258 6300

Sequencing Primers

Primer	Sequence
Forward	5'-gatccgcgagatcctgat-3'
Reverse	5'-cctattggcgttactatg-3'

GeneCopoeia, Inc.
9620 Medical Center Drive, #101
Rockville, MD 20850
Tel: 301-762-0888 Fax: 301-762-3888
Email: inquiry@genecopoeia.com Web: www.genecopoeia.com