


pEZX-GA01

6959bp, circular



MCS1: BstBI, BglII, EcoRI, EcoRV, HindIII
MCS2: NotI

Vector Features

Name	Position
Gluc	8 565
MCS2	566 573
SV40 polyA	589 710
CMV promoter	869 1378
SeAP	1442 2986
SV40 promoter	3227 3556
Puromycin resistance gene	3618 4217
SV40 polyA	4317 4438
Kanamycin resistance gene	4650 5459
pUC ori	5549 6137
TPS	6660 6884
MCS1	6885 6941

Sequencing Primers

Primer	Sequence
Forward	5'-gtcattctattctggggg-3'
Reverse	5'-ttgttctcggtgggcttggc-3'

GeneCopoeia, Inc.
9620 Medical Center Drive, #101
Rockville, MD 20850
Tel: 301-762-0888 Fax: 301-762-3888
Email: inquiry@genecopoeia.com Web: www.genecopoeia.com

GeneCopoeia products are for research use only

© 2016, GeneCopoeia, Inc