


pEZX-LvGA01

9732bp, circular



MCS1: AfeI, AgeI, PmeI, EcoRI, MfeI, BstBI
MCS2: XhoI, PspXI

Vector Features

Name	Position
SV40 promoter	1 368
SeAP	396 1940
IRES	1963 2536
Puromycin resistance gene	2569 3168
WPRE	3182 3770
3'LTR	3977 4210
pUC ori(complementary strand)	4730 5315
Ampicillin resistance gene(complementary strand)	5486 6346
SV40 polyA	6499 6633
5'LTR	6761 7172
HIV-1 psi pack sequence	7217 7342
RRE	7843 8076
cPPT	8568 8691
Synthesis polyA	8815 8868
MCS2	8869 8874
mGluc(complementary strand)	8875 9432
MCS1	9433 9470
TPS(complementary strand)	9471 9695

Sequencing Primers

Primer	Sequence
Forward	5'-ttgttctcgggtgggcttggc-3'
Reverse	5'-caccctaactgacacacattccac-3'

GeneCopoeia, Inc.
9620 Medical Center Drive, #101
Rockville, MD 20850
Tel: 301-762-0888 Fax: 301-762-3888
Email: inquiry@genecopoeia.com Web: www.genecopoeia.com

GeneCopoeia products are for research use only

© 2016, GeneCopoeia, Inc