

SPQ


Catalog Number	Packaging Size
C274	50 mg

Storage upon receipt: -20°C, protected from light

Introduction

The chloride indicator **SPQ** has absorption/emission maxima ~344/443 and detects Cl⁻ via diffusion-limited collisional quenching.

Specifications

Label:	SPQ	
Ex/Em:	344/443 nm	
Detection Method:	Fluorescent	
Solubility:	DMSO, DMF	
Molecular Formula	C ₁₃ H ₁₅ NO ₄ S	
Molecular Weight:	281.32	
CAS Number:	83907-40-8	
Storage Conditions:	-20°C, protect from light	
Shipping Condition:	Room Temperature	

Applications

Chloride indicator

References:

NO decreases thick ascending limb chloride absorption by reducing Na(+)-K(+)-2Cl(-) cotransporter activity.

Ortiz PA, Hong NJ, Garvin JL
Am J Physiol Renal Physiol (2001) 281:F819-F825

Endogenous nitric oxide inhibits chloride transport in the thick ascending limb.

Plato CF, Stoos BA, Wang D, Garvin JL
Am J Physiol (1999) 276:F159-F163

Characterisation of Chloride Transport at the Tonoplast of Higher Plants Using a Chloride-Sensitive Fluorescent Probe.

Pope AJ, Leigh RA
Planta (1990) 181:406-406