

CMV-SV40T (Puro) Lentifect™ Purified Lentiviral Particles
Cat No. LPP-SV40T-Lv105-100-C

Ready-to-use lentiviral particles for the transduction of a variety of mammalian cells including difficult-to-transfect, primary, stem and non-dividing cells.

Description

GeneCopoeia Lentifect™ Lentiviral Particles are produced from a standardized protocol using purified plasmid DNA and the proprietary reagents, EndoFectin™ Lenti (for transfection) and TiterBoost™ solution. The protocol uses a third generation self-inactivating packaging system meeting BioSafety Level 2 requirements.

The Lentifect particles include a CMV promoter for efficient expression of non-tagged, SV40 large T antigen protein in target cells and use a puromycin resistance marker for selection of stably transduced cells.

Contents and storage

Provided as 4 vials of 25 µl of purified SV40T lentiviral particles with titers of $\sim 1 \times 10^8$ TU/ml.

Lentifect particles are shipped on dry ice and must be stored at -80°C immediately upon receipt. Avoid repeated freeze-thaw cycles as this will reduce titers.

Quality control

The lentiviral expression construct was validated by full-length sequencing, restriction enzyme digestion and PCR-size validation using gene-specific and vector-specific primers. Product is confirmed free of bacteria, fungi and common Mycoplasma contamination.

Viral titer

The transduction unit (TU or IFU) of the lentiviral particles was estimated using the formula- $1\text{TU}=100$ copies of viral genomic RNA. The physical copy numbers of the viral genomic RNA was determined using qRT-PCR. The customer should test the transduction at MOI=0.3, 1, 3, 5, 10 for their specific cell lines in order to get the best transduction efficiency.

Overview of production

The SV40 Large T Antigen OmicsLink™ ORF lentiviral expression plasmid (GeneCopoeia Cat. No. EX-SV40T-Lv105) was constructed using GeneCopoeia proprietary RecJoin™ technology. This plasmid was co-transfected into 293Ta cells (GeneCopoeia Cat. No. LT008) with the Lenti-Pac HIV Packaging Mix (GeneCopoeia Cat. No. LT001-01). Lentivirus-containing supernatants were harvested 48 hours after transfection and stored at -80°C .

User manual

Please contact GeneCopoeia for a copy or download at:

http://genecopoeia.com/product/lentiviral/pdf/packaging_kit_manual.pdf

GeneCopoeia, Inc.
9620 Medical Center Drive, Suite 101
Rockville, Maryland 20850
Tel: 301-762-0888 Fax: 301-762-8333
Email: inquiry@genecopoeia.com
Web: www.genecopoeia.com

GeneCopoeia Products are for Research Use Only Copyright © 2019 GeneCopoeia Inc.
Trademarks: Lentifect™, Lenti-Pac™, OmicsLink™, EndoFectin™, TiterBoost™ (GeneCopoeia Inc.) LPPSV40T042619 www.genecopoeia.com