

Flag-copGFP-CD63 Lentifect™ Purified Lentiviral Particles • **Catalog #: LP786-025, LP786-100**

Lentiviral particles pseudotyped with VSV-G protein for transferring genes to a variety of mammalian cells including primary, stem and non-dividing cells.

Description

- Gene: Flag-copGFP-CD63
- Promoter: EF1a
- Tag: 3xFLAG
- Reporter: copGFP
- Resistance marker: Blasticidin resistance gene
- Additional note: N/A

Contents and storage

Provided as 1 vial of 25 µl or 4 vials of 25 µl of purified CD63 lentiviral particles with titers of $\sim 1 \times 10^8$ TU/ml.

Lentiviral particles are shipped on dry ice and must be stored at -80°C immediately upon receipt. Avoid repeated freeze-thaw cycles as this will reduce titers. Divide into useful aliquots if necessary and store at -80°C .

Viral titer

The transduction unit (TU or IFU) of the lentiviral particles was estimated using the formula- $1\text{TU}=30$ copies of viral genomic RNA. The physical copy numbers of the viral genomic RNA was determined using qRT-PCR. The customer should test the transduction at MOI=0.3, 1, 3, 5, 10 for their specific cell lines in order to get the best transduction efficiency.

Overview of production

Identity of lentiviral transfer vector: pEZ-Lv197. The lentiviral particles were generated by following a standardized protocol using highly purified plasmids and EndoFectin-Lenti™ and TiterBoost™ reagents. The lentiviral transfer vector was co-transfected into 293Ta cells (Cat #: LT008) with Lenti-Pac™ HIV packaging mix (Cat #: LT001). The lentiviral particles were purified and stored at -80°C in aliquots (purified particles).

Quality control

The lentiviral expression construct was validated by full-length sequencing, restriction enzyme digestion and PCR-size validation using gene-specific and vector-specific primers. Product is confirmed free of bacteria, fungi and common Mycoplasma contamination.

User manual

Please contact GeneCopoeia for a copy or download at:

<https://www.genecopoeia.com/wp-content/uploads/2018/03/Lentivirus-protocol-GeneCopoeia.pdf>

GeneCopoeia, Inc.

9620 Medical Center Drive, #101

Rockville, Maryland 20850

Tel: 301-762-0888 Fax: 301-762-8333

Email: inquiry@genecopoeia.com Web: www.genecopoeia.com

GeneCopoeia Products are for Research Use Only

Trademarks: EndoFectin™, Lenti-Pac™, TiterBoost™ (GeneCopoeia Inc.)

Copyright© 2023 GeneCopoeia Inc.

LP-20231212