

## Lentiviral particles for Firefly Luciferase

# Catalog #: SP005-400

Lentiviral particles pseudotyped with Truncated SARS-CoV-2 spike glycoprotein (D614G & N439K) for transferring genes to a variety of mammalian cells expressing human ACE2 receptor on the cell surface.

### Description

- Gene: Firefly Luciferase
- Promoter: CMV
- Tag: N/A
- Reporter : EGFP
- Resistance marker: Puromycin resistance gene
- Additional note: N/A

## Contents and storage

Provided as 50 µl/vial

Titer 1: > 5×10<sup>6</sup> TU/ml Titer 2: > 5×10<sup>8</sup> RLU/ml

#### For lot-specific datasheet and titers, please visit http://www.igenebio.com/tech/datasheet/index.php

Lentivirus particles are shipped on dry ice and must be stored at -80  $^{\circ}$ C immediately upon receipt. Avoid repeated freeze-thaw cycles as this will reduce titers. Divide into useful aliquots if necessary and store at -80  $^{\circ}$ C.

#### Viral titer

The titers of Transduction Unit/ml (TU/ml) and Relative Fluorescent Unit/ml (RFU/ml) of the spike protein-pseudotyped lentiviral particles were provided for each particle, based on the transduction efficiency using GeneCopoeia's ACE2 and TMPRSS2/HEK293T stable cell line (Cat.# SL222).

The luciferase titer (RLU/ml) was based on the viral particle's firefly luciferase signal after transducing the ACE2 and TMPRSS2/HEK293T cells, and was determined using GeneCopoeia's Luc-Pair<sup>™</sup> Firefly Luciferase HS Assay Kit (Cat.# LF007). The eGFP transduction titer (TU/ml) was determined by the FACS sorting of the eGFP expressing ACE2 and TMPRSS2 /HEK293T cells after pseudotyped viral transduction.

The customer should test the transduction at MOI=0.3, 1, 3, 5, 10 for their specific cell lines in order to get the best transduction efficiency.

#### **Overview of production**

Identity of lentiviral transfer vector: pEZ-Lv201

The lentiviral particles were generated by following a standardized protocol using highly purified plasmids and

EndoFectin-Lenti<sup>™</sup> and TiterBoost<sup>™</sup> reagents.

The lentiviral transfer vector was co-transfected into 293Ta cells (Cat #: LT008) with Lenti-Pac SARS-CoV-2 Truncated Spike protein (D614G & N439K) packaging mix(Cat#: LT033). The lentivirus particles were purified and stored at -80 °C in aliquots (purified particles).

## **Quality control**

The lentiviral expression construct was validated by full-length sequencing, restriction enzyme digestion and PCR-size validation using gene-specific and vector-specific primers. Product is confirmed free of bacteria, fungi and common Mycoplasma contamination.

#### User manual

Please contact GeneCopoeia for a copy or download at: Lenti-Pac<sup>™</sup> SARS-CoV-2 Spike protein-pseudotyped Lentivirus Packaging Kit

GeneCopoeia, Inc. 9620 Medical Center Drive, #101 Rockville, Maryland 20850 Tel: 301-762-0888 Fax: 301-762-8333 Email: <u>inquiry@genecopoeia.com</u> Web: <u>www.genecopoeia.com</u>

www.genecopoeia.com