


# OmicLink™ Expression Clone (CMV Promoter)



**Vector: pEZ-M98 7982 bps DNA Circular**

**Vector Features :**

attR2: bases 7866-8  
 monomeric eGFP: bases 28-744  
 BGH polyadenylation region: bases 964-1191  
 f1 origin: bases 1237-1665  
 SV40 early promoter and origin : bases 1692-2001  
 Neomycin resistance ORF: bases 2075-2869  
 SV40 early polyadenylation region: bases 3046-5175  
 pUC origin: bases 3558-4267  
 Ampicillin resistance ORF: bases 4376-5236 (complementary strand)  
 CMV promoter: bases 5604-6191  
 T7 promoter: bases 6235-6254  
 attR1: bases 6286-6410  
 CmR: bases 6519-7178  
 ccdB: bases 7520-7825