


EZShuttle™ Recombination Cloning Expression Vector pEZ-B01

7115 bp, circular



Vector Features

Name	Start	End
T7 terminator	490	537
f1 origin	574	1029
Ampicillin resistance gene (complementary strand)	1144	2004
AmpR promoter (complementary strand)	2005	2097
pUC ori	2178	2766
<i>lacI</i> (complementary strand)	3959	5041
T7 promoter	5432	5450
lac operator	5451	5475
6xHis	5523	5540
<i>attR1</i>	5547	5671
<i>ccdB</i>	6769	7074
<i>attR2</i>	7115	26

Sequencing Primers

Primer	Sequence
Forward	5'- taatacgactcactataggg -3'
Reverse	5'- ttcacttctgagttcggcatg -3'

GeneCopoeia, Inc.
9620 Medical Center Drive, #101
Rockville, MD 20850
Tel: 301-762-0888 Fax: 301-762-3888
Email: inquiry@genecopoeia.com Web: www.genecopoeia.com

GeneCopoeia products are for research use only

© 2015, GeneCopoeia, Inc