


## EZShuttle™ Recombination Cloning Expression Vector pEZ-M06

7549bp, circular



### Vector Features

Name	Position
T7 promoter(complementary strand)	42 60
SV40 poly A	358 625
f1 ori	751 1005
SV40 promoter	1366 1675
Neomycin resistance gene	1733 2527
Synthesis poly A	2594 2642
bla promoter	2953 3051
Ampicillin resistance gene	3052 3912
pUC origin	4058 4731
CMV promoter	5132 5725
3xHA tag	5751 5837
attR1	5853 5977
lac UV5 promoter	6002 6032
CmR	6086 6745
ccdB	7087 7392
attR2	7433 8

### Sequencing Primers

Primer	Sequence
Forward	5'- cagcctccggactctagc -3'
Reverse	5'- taatagactcactataggg-3'

GeneCopoeia, Inc.  
9620 Medical Center Drive, #101  
Rockville, MD 20850  
Tel: 301-762-0888 Fax: 301-762-3888  
Email: [inquiry@genecopoeia.com](mailto:inquiry@genecopoeia.com) Web: [www.genecopoeia.com](http://www.genecopoeia.com)

GeneCopoeia products are for research use only

© 2015, GeneCopoeia, Inc