


pEZX-GA02

8577bp, circular



MCS2: AsiSI, EcoRI, BsiWI, PmeI, XhoI, SpeI

Vector Features

Name	Position
SV40 polyA	51 242
CMV promoter	292 819
SeAP	912 2456
bGH pA	2578 2805
f1 ori	2851 3279
SV40 promoter	3306 3615
Neomycin resistance gene	3689 4483
SV40 polyA	4660 4789
pUC ori (complementary strand)	5172 5845
Ampicillin resistance gene (complementary strand)	5990 6850
SV40 promoter	7282 7701
chimeric intron	7764 7896
Gluc	7977 8534
MCS2	8535 8577

Sequencing Primers

Primer	Sequence
Forward	5'-atctgcctgtcccacatcaa-3'
Reverse	5'-CCTATTGGCGTTACTATG-3'

GeneCopoeia, Inc.
9620 Medical Center Drive, #101
Rockville, MD 20850
Tel: 301-762-0888 Fax: 301-762-3888
Email: inquiry@genecopoeia.com Web: www.genecopoeia.com