


EZShuttle™ Recombination Cloning Expression Vector

pEZ-Lv233.1

10846 bp, circular



Vector Features

Name	Start	END
SV40 promoter	17	330
EGFP	410	720
IRES	1146	1719
PuroR	1752	2351
WPRE	2365	2593
3'LTR (ΔU3)	3160	3393
pUC ori	3913	4498
AmpR	4469	5529
5' LTR	6173	6353
HIV-1 RRE	7026	7259
cPPT/CTS	7751	7866
EF-1α promoter	7971	9148
attR1	9155	9275
CmR	9338	10047
ccdB	10377	10682
attR2	10723	10846

Sequencing Primers

Primer	Sequence
Forward	5'-GAGGGGTTTTATGCGATG-3'
Reverse	5'-CTGGAATAGCTCAGAGGC-3'

GeneCopoeia, Inc.
9620 Medical Center Drive, #101
Rockville, MD 20850
Tel: 301-762-0888 Fax: 301-762-3888
Email: inquiry@genecopoeia.com Web: www.genecopoeia.com

GeneCopoeia products are for research use only

© 2020, GeneCopoeia, Inc