

EndoFectin™ 转染试剂

高效、低毒、适用于多种细胞系

买一赠一

活动时间：2025 年 3 月 1 日 ~ 6 月 30 日



扫码查看产品详情



- 转染效率高，适用于多种细胞系的转染
- 细胞毒性低，保证高细胞活性
- 转染前无需去除细胞培养液或血清
- 转染后无需清洗细胞或更换培养基
- 可用于基因过表达、基因敲减/敲除、微孔板的高通量转染操作




图 . Max 转染试剂简便操作流程

明星产品

产品名称	适用性	货号	规格	目录价	促销价
EndoFectin™ Max Hot!	广泛适用于多种细胞	EF013	1 mL	¥ 2000	买一赠一
		EF014	3 mL	¥ 5000	买一赠一

* 凡同一订单购买 EndoFectin™ Max，加购 miTarget™ miRNA 靶标 (3' UTR) 克隆或 GLuc-ON™ 启动子报告克隆，则该克隆享 9 折优惠，且该 EndoFectin™ Max 在 7 折促销基础上再享额外 8 折。

EndoFectin™ 系列除了明星产品 Max，更有针对 293T、Cas9 蛋白转染的专业转染试剂！

产品名称	适用性	货号	规格	目录价	促销价
EndoFectin™ Lenti	转染 293T 细胞，可收获更高滴度的慢病毒	EF001	1 mL	¥ 1700	买一赠一
		EF002	3 mL	¥ 4000	买一赠一
CRISPR-Fectin™	CRISPR-Cas9 核蛋白专用的转染试剂，可与 GeneHero™ Cas9 蛋白配合使用	EF015	1mL	¥ 2000	买一赠一
		EF016	3mL	¥ 5000	买一赠一
EndoFectin™ Expi293	专用于分泌蛋白表达 专用于胞内蛋白表达	EF017	1L 转染用量	¥ 3300	买一赠一
		EF019	1L 转染用量	¥ 3300	买一赠一
EndoFectin™ RNAi	专为将 siRNA 和 miRNA 递送至哺乳动物细胞而设计	EF021	1 mL	¥ 2000	买一赠一
		EF022	3 mL	¥ 5000	买一赠一

客户发表文章举例：

- ◆ The human mitochondrial translation factor TACO1 alleviates mitoribosome stalling at polyproline stretches. 2024. *Nucleic Acids Research*. **EndoFectin™ Max.**
- ◆ Mitoribosome structure with cofactors and modifications reveals mechanism of ligand binding and interactions with L1 stalk. 2024. *Nature Communications*. **EndoFectin™ Max.**
- ◆ Hyperglycemia-Enhanced Neutrophil Extracellular Traps Drive Mucosal Immunopathology at the Oral Barrier. 2024. *Advanced Science*. **EndoFectin™ RNAi.**