

Plate	Position	Catalog No. of Primer	Accession No. of Gene	Symbol
QG086-01	A01	HQP008832	NM_181054	HIF1A
QG086-01	A02	HQP009026	NM_000194	HPRT1
QG086-01	A03	HQP018485	NM_000551	VHL
QG086-01	A04	HQP021496	NM_003807	TNFSF14
QG086-01	A05	HQP018148	NM_001065	TNFRSF1A
QG086-01	A06	HQP018141	NM_000594	TNF
QG086-01	A07	HQP018116	NM_138554	TLR4
QG086-01	A08	HQP018114	NM_003264	TLR2
QG086-01	A09	HQP022584	NM_004787	SLIT2
QG086-01	A10	HQP016672	NM_006923	SDF2
QG086-01	A11	HQP022410	NM_004757	SCYE1
QG086-01	A12	HQP011807	NM_003998	NFKB1
QG086-01	B01	HQP011603	NM_002468	MYD88
QG086-01	B02	HQP011258	NM_002423	MMP7
QG086-01	B03	HQP009679	NM_000634	IL8RA
QG086-01	B04	HQP009678	NM_000584	IL8
QG086-01	B05	HQP009662	NM_000589	IL4
QG086-01	B06	HQP009640	NM_000575	IL1A
QG086-01	B07	HQP009718	NM_001562	IL18
QG086-01	B08	HQP009697	NM_002188	IL13
QG086-01	B09	HQP008831	NM_001530	HIF1A
QG086-01	B10	HQP007585	NM_032554	GPR81
QG086-01	B11	HQP008143	NM_005299	GPR31
QG086-01	B12	HQP020039	NM_000557	GDF5
QG086-01	C01	HQP004538	NM_001953	ECGF1
QG086-01	C02	HQP054015	NM_014376	CYFIP2
QG086-01	C03	HQP000808	NM_006564	CXCR6
QG086-01	C04	HQP018803	NM_003467	CXCR4

QG086-01	C05	HQP007900	NM_001504	CXCR3
QG086-01	C06	HQP011220	NM_002416	CXCL9
QG086-01	C07	HQP016648	NM_002993	CXCL6
QG086-01	C08	HQP016650	NM_002994	CXCL5
QG086-01	C09	HQP008459	NM_002090	CXCL3
QG086-01	C10	HQP008458	NM_002089	CXCL2
QG086-01	C11	HQP000672	NM_006419	CXCL13
QG086-01	C12	HQP016669	NM_000609	CXCL12
QG086-01	D01	HQP016649	NM_005409	CXCL11
QG086-01	D02	HQP009746	NM_001565	CXCL10
QG086-01	D03	HQP008456	NM_001511	CXCL1
QG086-01	D04	HQP003691	NM_001337	CX3CR1
QG086-01	D05	HQP016652	NM_002996	CX3CL1
QG086-01	D06	HQP003173	NM_000759	CSF3
QG086-01	D07	HQP002242	NM_004072	CMKLR1
QG086-01	D08	HQP003295	NM_178818	CMTM4
QG086-01	D09	HQP054031	NM_144601	CMTM3
QG086-01	D10	HQP003297	NM_144673	CMTM2
QG086-01	D11	HQP001619	NM_181269	CMTM1
QG086-01	D12	HQP012497	NM_181641	CKLF
QG086-01	E01	HQP021922	NM_003965	CCRL2
QG086-01	E02	HQP012851	NM_016557	CCRL1
QG086-01	E03	HQP002222	NM_005201	CCR8
QG086-01	E04	HQP002217	NM_001838	CCR7
QG086-01	E05	HQP002212	NM_004367	CCR6
QG086-01	E06	HQP002210	NM_000579	CCR5
QG086-01	E07	HQP002209	NM_005508	CCR4
QG086-01	E08	HQP002201	NM_000648	CCR2
QG086-01	E09	HQP007833	NM_016602	CCR10
QG086-01	E10	HQP002198	NM_001295	CCR1

QG086-01	E11	HQP016628	NM_005623	CCL8
QG086-01	E12	HQP016627	NM_006273	CCL7
QG086-01	F01	HQP016626	NM_002985	CCL5
QG086-01	F02	HQP016625	NM_002984	CCL4
QG086-01	F03	HQP016622	NM_002983	CCL3
QG086-01	F04	HQP016621	NM_002982	CCL2
QG086-01	F05	HQP016638	NM_006274	CCL19
QG086-01	F06	HQP016637	NM_002988	CCL18
QG086-01	F07	HQP016636	NM_002987	CCL17
QG086-01	F08	HQP016635	NM_004590	CCL16
QG086-01	F09	HQP054027	NM_032965	CCL15
QG086-01	F10	HQP016630	NM_005408	CCL13
QG086-01	F11	HQP016629	NM_002986	CCL11
QG086-01	F12	HQP016620	NM_002981	CCL1
QG086-01	G01	HQP002229	NM_001296	CCBP2
QG086-01	G02	HQP018326	NM_001736	C5AR1
QG086-01	G03	HQP018304	NM_001735	C5
QG086-01	G04	HQP016959	NM_001716	BLR1
QG086-01	G05	HQP004535	NM_005161	AGTRL1
QG086-01	G06	HQP002207	NM_001837	CCR3
QG086-01	G07	HQP001625	NM_181296	CMTM1
QG086-01	G08	HQP016545	NM_001709	BDNF
QG086-01	G09	HQP009714	NM_004513	IL16
QG086-01	G10	HQP054072	NM_004530	MMP2
QG086-01	G11	HQP017986	NM_004610	TCP10
QG086-01	G12	HQP013418	NM_018643	TREM1
QG086-01	H01	HGDC		
QG086-01	H02	HGDC		
QG086-01	H03	HQP006940	NM_002046	GAPDH
QG086-01	H04	HQP016381	NM_001101	ACTB

QG086-01	H05	HQP015171	NM_004048	B2M
QG086-01	H06	HQP006171	NM_012423	RPL13A
QG086-01	H07	HQP009026	NM_000194	HPRT1
QG086-01	H08	HQP054253	NR_003286	RN18S1
QG086-01	H09	RT		
QG086-01	H10	RT		
QG086-01	H11	PCR		
QG086-01	H12	PCR		